

M A X I C O M F O R T TM S E R I E S



**Pocketed Coil
Spring Seat
Standard!***



Bourbon



Geranium



Coffee Bean*



Bridle*

*** upcharge applies**



**Pub
PR713
shown in Bridle**

**Pub
PR713**

The Pub chair is a traditional style two pillow waterfall back recliner featuring decorative brass nail head trim on both the arms and wings. The supple waterfall, zipper back design allows for customizing comfort. This dual motor chair is part of the MaxiComfort series. Recline in style!



Pub Chair

	Pub
Model Number	PR713
Size	Medium
Weight Capacity	375 Lb.
Positioning	MaxiComfort
Chaise Seat	YES
Pocketed Coil Seat	YES
Back Type	Pillow
Removable Back Option	NO
Head Rest Cover	NO
Arm Rest Covers	NO
Battery Back Up	YES
Overall Width	35"
Overall Height	45"
Width Between Arms	21"
Floor to Top of Seat	20"
Seat Depth	18"
Seat to Top of Back	26.5"
Seat to Top of Arm	4.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	23"
Distance Required from Wall Reclined	25"
Extension reclined from Back of Seat to End of Foot Rest	43.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Option	NO
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option	NO
Extra Pocket Option	Standard
Heat & Massage Option	YES
Power Pillow Option	NO
Power Adjustable Lumbar Option	NO
AutoFlex Option	NO
Foot Rest Extension Option	NO
Straight Lift Option	NO
USB Charging Port	NO
Place of Origin	USA



Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
UL 73 and
CSA C22.2 No. 68

