

The Essential Collection from Pride® Power Lift Recliners provides exceptional value and outstanding comfort. Several fabric options will match any décor while delivering rest and worry-free operation.



Easy-to-use hand control

ENERGY EFFICIENT WITH BUILT-IN BATTERY BACKUP



Scan here for more information on the **Essential Collection**®



FOLLOW US ON SOCIAL MEDIA

FEATURES & BENEFITS


- Easy-to-use hand control adjusts your chair to find your perfect position
- Large dual pockets offer ample storage space that is easily accessible
- Smooth and durable lift and scissor mechanism for long-lasting and quiet performance
- Removable back lets you adjust chair padding as needed
- Integrated battery backup ensures chair is operational when the power goes out
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Lifetime warranty on steel lift frame
- Fire retardant (California specifications; bulletin 117-2013)

OPTIONS*

- Cleaning kit
- Non-skid leg levelers

* Additional charges apply for all options. Extra lead time may apply for some options.

MODEL	LC-250	LC-105
POSITION	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'2" - 5'10"
SIZE	Medium	Medium
WEIGHT CAPACITY	375 lbs.	325 lbs.
SEAT-TO-FLOOR	20"	18"
SEAT DEPTH	20"	19"
SEAT WIDTH	22"	20"
TOP OF BACK TO SEAT	26"	28"
FULLY PADDED CHAISE	Yes	Yes
SEATING TYPE	Total Comfort	Sinuous Spring
DISTANCE FROM WALL	20"	17"
OVERALL WIDTH	34"	33"
CHAIR WEIGHT	118 lbs.	96 lbs.
WARRANTY	2-Year Limited	1-Year Limited



L-250 FABRICS

LC-105 FABRICS

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]